


**AVATAR<sup>NX</sup>**  
UNCEMENTED HIP SYSTEM

Designed to accurately restore joint biomechanics & address the demand for bone conservation




Stem Size	1	2	3	4	5	6 <sup>#</sup>	7 <sup>#</sup>	8 <sup>#</sup>
Length	130	140	146	150	154	161	166	170
Head Offset	39	39	41	41	42	42	43	44
Neck Length	38.5							
Neck Angle	135°							

Stem length, head offset & neck length in mm

\* Note: head offset mentioned is for sizes (mm): 28, 32, 36 (0)  
# Will be available on request

Designed for cemented  
implantation



## Bipolar head

### UHMWPE liners

Helps to reduce wear

### Pre-assembled monoblock construct

Prevents micro-motion & allows intra-operative ease of use

### Highly polished SS outer shell

Designed to reduce friction

### Easy to use

### efficient instruments

Designed to enhance efficiency

### Self - centering design

## Femoral head



### Multiple neck length head options

To optimally restore joint biomechanics intra-op

### Internal head sizes (mm)

For OD: 37 - 43 = 22.2 (-2.5, 0, +2.5)

For OD: 45 - 53 = 28 (-4, 0, +4)

(2mm Increments In OD)

\* Note: all 12/14 heads available in the bioradmedisys portfolio are compatible with the Avatar NX uncemented and cemented stems.